

KRAS/BRAF Mutation Analysis

ONCOLOGY

KRAS/BRAF Mutations and Cancer

The *KRAS* gene encodes a small GTPase that plays a key role in transducing signals from the epidermal growth factor receptor (EGFR) to downstream effectors. *KRAS* mutations have been commonly found in several types of human malignancies, such as metastatic colorectal cancer (mCRC), lung adenocarcinoma and thyroid cancer. Activating mutations are commonly found in exon 2 (codon 12/13), exon 3 (codon 61), and exon 4 (codon 117/146).

Several studies have demonstrated that tumors carrying any of these mutant forms of the *KRAS* gene are less likely to respond to anti-EGFR antibody therapy. The American Society of Clinical Oncology (ASCO) recently released its first Provisional Clinical Opinion (PCO) suggesting that all patients to be administered anti-EGFR monoclonal antibody and tyrosine kinase inhibitor therapies (e.g. cetuximab, panitumumab, gefitinib and erlotinib) should be screened for *KRAS* mutations.

Recent studies have also shown that not all mCRC patients with wild-type *KRAS* tumors respond to anti-EGFR therapy. This suggests that additional genes and/or pathways may be involved in the mechanism of resistance to these drugs. Mutations in exon 4 of *KRAS*, as well as exon 15 of *BRAF*, another downstream effector of the EGF-activated pathway, have been identified in up to 15% of mCRC tumors. Patients with these activating mutations have decreased progression-free (PFS) and overall (OS) survival when treated with EGFR antagonists.

These findings strongly suggest that screening for both *KRAS* and *BRAF* mutations is necessary to more accurately identify patients who will not respond to anti-EGFR therapy.

Testing Procedure and Analysis

EntroGen's *KRAS/BRAF* mutation panel is a polymerase chain reaction (PCR)-based assay and uses allele-specific primers to identify the presence of mutations in *KRAS* codons 12, 13, 61, 117 and 146 as well as *BRAF* V600E (also available separately). The testing procedure involves three (3) simple steps:

- Isolation of DNA from tumor biopsies, paraffin-embedded sections (FFPE), or fresh frozen tumors
- Amplification of regions of the *KRAS* and *BRAF* genes using allele-specific primers
- Detection of amplification product on a Real-Time PCR instrument.

This test can be completed in approximately 2 hours from isolation of DNA to test result.

Equipment and Materials

EntroGen's *KRAS/BRAF* mutation panel requires a real-time PCR instrument capable of detecting FAM and VIC fluorescent probes.

This test includes reagents required for the PCR amplification/detection, as well as validated reaction controls. Columns and reagents for DNA isolation are not included.

Intended Use

EntroGen's *KRAS/BRAF* mutation analysis kits are available for Research Use Only (RUO) in the United States and for In Vitro Diagnostic use (CE-IVD) in Europe.

Ordering Information

Product Name	Cat. No.
K-Ras Mutation Analysis Kit for Real-Time PCR (exons 2, 3 & 4)	KRAS-RT50
B-Raf Mutation Analysis Kit for Real-Time PCR (exon 15)	BRAF-RT50
K-Ras/B-Raf Mutation Analysis Kit for Real-Time PCR (exons 2, 3 & 4 of <i>KRAS</i> and exon 15 of <i>BRAF</i>)	KRBR-RT50